





Cerebral Open Flow Microperfusion (cOFM) for In-Vivo Cerebral Fluid Sampling -Comparison of cOFM and Microdialysis

October 14th, 2021



#### Who We Are

# Today's Speaker



#### Joanna Hummer, PhD

Deputy Head of Biomedical Tissue Monitoring Joanneum Research HEALTH

joanna.hummer@joanneum.at



#### Florie Le Prieult, PhD

Lab Head - In vivo Neuroscience Pharmacokinetics AbbVie

florie.leprieult@abbvie.com



#### Who We Are

#### Boutique CRO for research projects & drug development programs



# We support clinical and preclinical activities by providing:

- (Pre)clinical PK/PD/BE studies at the target tissue level
- In-vitro release testing (IVRT)
- Customized bioanalyses (PK, PD) GLP certified lab
- Metabolomics
- Data management
- Biostatistics
- Medical writing

#### Who We Are

# Open Flow Microperfusion (cOFM)










- ✓ provides unfiltered interstitial fluid
- / PK/PD at target tissue
- ✓ CE-certified for clinical use (aOFM/dOFM)






# Biomedical Tissue Monitoring in the Brain



#### Biomedical Tissue Monitoring

#### How to Measure Substances in the Brain In-Vivo?



Credit: https://www.researchgate.net/figure/Brain-vasculature-separates-the-circulating-blood-from-the-CNS-tissue-The-vessels-are\_fig1\_320026611

Credit: JOANNEUM RESEARCH




#### Question

# I am working with ...

- MD
- a/dOFM
- cOFM
- others






# Cerebral Open Flow Microperfusion




# Concept




cOFM is a unique tool to investigate pharmacokinetics (PK) and pharmacodynamics (PD) in the brain with an intact BBB.



# Implantation and Sampling



- 1 Guide Ø 0.5 mm
- 2 Healing Dummy
- 3 Sampling Insert
- 4 Locking Wedge



Implantation and sampling are 14 days apart

⇒ intact Blood Brain Barrier




#### Intact Blood Brain Barrier

#### **Evans Blue BBB permeability**

# BBB permeability for EB after cOFM implantation Control group without EB and a control group without EB


#### BBB marker – Sodium fluorescein (Naf)








# Time-Resolved Substance Concentration Monitoring






#### Minimal Tissue Reaction

#### **Quantification of Astrocytes**



#### **Quantification of Macrophages**



Distance to cOFM probe Distance to cOFM probe

- ✓ cOFM probe material elicits minimal tissue reaction.
- √ No glial scar
- ✓ 15 days and 30 days after cOFM implantation





# Different Brain Regions


#### **ISF** sampling in

- striatum
- cortex
- hippocampus
- hypothalamus
- ..and others

#### **CSF** sampling in

lateral ventricle

third ventricle





# Monitoring of Substances without Limitation

#### Lipophilic substances

Amitriptyline (logP: 4.9), fluoxetine (logP: 4.6), ...

#### Big substances

- Antibodies (trastuzumab, anti-BACE1, ...)
- Proteins (amyloid ß, tau, albumin, IGGs, …)
- Nanoparticles (doxorubicin, liposomes, ...)

#### Hormones

Leptin (16 kDa), GLP-1 (3297 Da)

#### Biomarkers

Cytokines, eicosanoids, growth factors, glucose, ...





#### Question

# Proteins or compounds that I am interested in measuring are:

- a) Neurotransmitters
- b) Neuropeptides
- c) Hormones
- d) Antibodies
- e) Drugs
- f) Nanoparticles
- g) Proteins
- h) Growth factors
- h) Cytokines
- J) others



# **Applications**

#### **BBB**

Monitoring of BBB function and permeability changes

#### **Pharmacokinetics**

Monitoring of substance transport across the intact BBB & time-resolved concentration profile in the brain

#### **Pharmacodynamics**

Monitoring of biomarkers in the brain

#### **Established:**

Rodents (mice, rats)

#### **Ongoing Development:**


- Pigs
- Non-human primates
- Humans





# Amitriptyline Sampling with cOFM and Microdialysis

Small (277 Da); lipophilic (logP 4.92); highly protein bound (95%); dosing: 25 mg/kg intraperitoneal



Analyte	Sampling Method	AUC (ng*h/mL)
AMI	cOFM	33.5
	MD	16.5
HYA	cOFM	9.4
	MD	n.a.
ANO	cOFM	13.9
	MD	n.a.
NOR	cOFM	n.a.
	MD	n.a.







#### Neurodegenerative Diseases

Collecting antibodies and large molecule biomarkers in mouse interstitial brain fluid: a comparison of microdialysis and cerebral open flow microperfusion. Le Prieult et al. (2021), mAbs, 13:1, 1918819, DOI: 10.1080/19420862.2021.1918819

#### Obesity/Diabetes

Time-resolved hypothalamic open flow micro-perfusion reveals normal leptin transport across the blood–brain barrier in leptin resistant mice, Maximilian Kleinert (2018). Molecular Metabolism, Volume 13, Pages 77-82. DOI: 10.1016/j.molmet.2018.04.008

#### Inflammation/BBB integrity

Assessment of blood-brain barrier function and the neuroinflammatory response in the rat brain by using cerebral open flow microperfusion (cOFM). Ghosh et al., PLoS One. 2014 May 22;9(5). DOI: 10.1371/journal.pone.0098143.

#### Nanoparticles

Enhanced doxorubicin delivery to the brain administered through glutathione PEGylated liposomal doxorubicin (2B3-101) as compared with generic Caelyx,(®) / Doxil (®) -- cerebral open flow microperfusion pilot study. Birngruber et al. (2014) J Pharm Sci. 2014 Jul;103(7):1945-1948. DOI: 10.1002/jps.23994

#### Glioblastoma (work in progress)



#### question

# I am interested in learning about

- a) cOFM & Neurodegeneration
- b) cOFM & Obesity/Diabetes
- c) cOFM & Glioblastoma
- d) cOFM & Nanoparticles
- E) cOFM & BBB integrity monitoring
- F) cOFM & Absolut quantification



### cOFM service @ JR HEALTH

- COFM setup is combined with Research PRODUCTS

  Raturn® and Culex®
- Simultaneous dosing & sampling (cISF, CSF & blood)
- Awake, freely moving animals









cOFM samples cerebral fluids with an intact blood brain barrier.



cOFM is membrane-free and samples all substances from the interstitial fluid independent of size, lipophilicity, or protein-binding.



cOFM enables continuous sampling for up to several days and a long implantation period for up to several weeks.



cOFM studies provide data for drug development as well as unique insights into brain metabolism and signaling.

#### Contact us

<u>joanna.hummer@joanneum.at</u> <u>cofm-services@joanneum.at</u>

#### Visit our websites

openflowmicroperfusion.com croservices.joanneum.at

#### Follow us



lcons in this presentation are made by Smashicons and Freepick from www.flaticon.com & edited by JOANNFUM RESEARCH HEALTH